

## Op-d-op II visor shields



- An advanced face shield with all the advantages of the original Op-d-op, with the added ability to *resize the visor at the click of a button!*
- Safety with comfort, style and convenience.
- Provides a perfect fit and eliminates uncomfortable pressure points.
- Fog resistant with increased ventilation.
- Change the crystal-clear face-shield in seconds.
- Available in 8 colours: Cloud white, Soft pink, Blue, Periwinkle, Pink, Violet, Green and Multi-colour.
- **SAVE £5 from our normal Visor pack price during MAY 2005.**

**Visor** Plus 3 replacement standard shields and 1 mini shield: **£19.00\***


**Face shields** Pack of 3 standard or 3 mini shields: **£10.00**

Pack of 12 standard or 12 mini shields: **£28.00**

\*Normal visor price £24.00. All prices subject to carriage charge and VAT

## Re-usable Bipolar & Diathermy forceps

We supply re-usable Bipolar and Diathermy forceps for all specialities. We also supply leads for both types of forceps. Contact us at the number below for a full colour brochure



DIA-021 Diathermy Forceps

BP-200-ISF Bipolar Forceps

## Surgical instrument cleaning brushes

Our comprehensive range of cleaning equipment is designed for use with all kinds of surgical instruments including cannulae, tubes, laparoscopic instruments, orthopaedic reamers and the Laryngeal Mask.

General instrument cleaning brush  
7.25" with stainless steel bristles for use on rasps, files reamers, burrs, needle holders, jaws, vascular serrations and stubborn stains.

General cleaning brush  
9" utility brush for use on sinks, autoclaves, ultrasonic cleaners, bowls, basins, I.V. poles, commodes, wheelchairs etc.



General instrument cleaning brushes  
single and double-ended, 7" and 7.25" with nylon bristles for use on scissors, hemostats, needle holders, dressing and tissue forceps.

Upright instrument cleaning brush  
5.25" with fine stainless steel bristles for use on rasps, files reamers, burrs, needle holders, jaws, vascular serrations and stubborn stains.

Upright instrument cleaning brush  
5.25" with nylon bristles for use on scissors, hemostats, needle holders, dressing and tissue forceps.

Examples from the wide range of cleaning brushes available



## High quality fibre optic light cables

**Prices for a complete cable, light source and instrument adaptor assembly start at just £225** *plus carriage and VAT*

Universal fibre optic cables with adaptors for all light sources and instruments. Up to 35 possible combinations.

- 6', 7.5' and 9' lengths
- 1,2,3,4,5 and 7mm diameters
- Special lengths available
- Colour range: Grey, Pink, Purple, Blue, Black, Red, Green and Gold
- All light source and instrument ends
- Autoclavable



## Colour coded locking tags

Shown here actual size

A practical and effective marking system of brightly coloured tags that quickly identify or mark instruments, carts, trays, cabinets and other equipment.

**Features:** Available in 12 distinctive colours • Non-sterile  
 Autoclavable • Tamper resistant • Latex Free  
 Durable and economical • Quick and easy to attach  
 Versatile - adapt to your own specific coding needs • Brightly coloured tags replace hard to see safety pins and sutures

**Click!**

Just £18.00 per 100 tags

## Bipolar and Monopolar Leads

Manufactured by **BOWA Electronic GmbH**, who use the technology they have developed to make the electrical cabling for BMW and Volkswagen cars. These cables, in their distinctive grey colouring, are **SEALED UNITS** and so avoid fluid ingress which is the chief cause of cable failure and deterioration.



**For more information** on the products featured in this Bulletin or for help with ordering, please call

**0845 310 1103**

All calls charged at local rates

## SIGH

Surgical Instrument Group Holdings Limited 89a Gloucester Road Croydon CR0 2DN  
 Telephone 020 8683 1103 Lo Call 0845 310 1103 Fax 020 8683 1105 Email [info@sightltd.com](mailto:info@sightltd.com)



[www.sightltd.com](http://www.sightltd.com)