

SIGH

BULLETIN

For Theatre and Sterile Services Managers

Surgical instrument cleaning brushes

Designed for use with all kinds of surgical instruments...

Our comprehensive range of cleaning equipment is designed for use with all kinds of surgical instruments including cannulae, tubes, laparoscopic instruments, orthopaedic reamers and the Laryngeal Mask.

Examples from the wide range of cleaning brushes available
Call 020 8683 1103 for more information

General instrument cleaning brush
7.25" with stainless steel bristles for use on rasps, files reamers, burrs, needle holders, jaws, vascular serrations and stubborn stains.

General cleaning brush
9" utility brush for use on sinks, autoclaves, ultrasonic cleaners, bowls, basins, I.V. poles, commodes, wheelchairs etc.



General instrument cleaning brushes
single and double-ended, 7" and 7.25" with nylon bristles for use on scissors, hemostats, needle holders, dressing and tissue forceps.


Upright instrument cleaning brush
5.25" with fine stainless steel bristles for use on rasps, files reamers, burrs, needle holders, jaws, vascular serrations and stubborn stains.

Upright instrument cleaning brush
5.25" with nylon bristle for use on scissors, hemostats, needle holders, dressing and tissue forceps.

Diathermy Quivers

Detailed leaflet available on request

- 95-010-10 Julian Taylor 20cm round base with drainage holes
- 95-010-20 20cm X 6cm flat base
- 95-010-30 Laparoscopic Quiver 38cm X 6cm



A complete fibre optic light cable, light source & instrument adaptor assembly for just £250

Universal fibre optic cables with adaptors for all light sources and instruments. Up to 35 possible combinations.

- 6', 7.5' and 9' lengths
- 1,2,3,4,5 and 7mm diameters

- Special lengths available
- Colour range: Grey, Pink, Purple, Blue, Black, Red, Green and Gold
- All light source and instrument ends
- Autoclavable



For more information, prices or to order any of these items - telephone 0845 310 1103



Shown here actual size

A practical and effective marking system of brightly coloured tags that quickly identify or mark instruments, carts, trays, cabinets and other equipment.

Features: Available in 12 distinctive colours • Non-sterile
 Autoclavable • Tamper resistant • Latex Free
 Durable and economical • Quick and easy to attach
 Versatile - adapt to your own specific coding needs • Brightly coloured tags replace hard to see safety pins and sutures

Click!

Just £18.00 per 100 tags



Re-usable Bipolar & Diathermy Forceps

We supply re-usable Bipolar and Diathermy forceps for all specialities. We also supply leads for both types of forceps. Contact us at the number below for a full colour brochure



DIA-021 Diathermy Forceps BP-200-ISF Bipolar Forceps

Bipolar and Monopolar Leads

Manufactured by **BOWA Electronic GmbH**, who use the technology they have developed to make the electrical cabling for BMW and Volkswagen cars. These cables, in their distinctive grey colouring, are **SEALED UNITS** and so avoid fluid ingress which is the chief cause of cable failure and deterioration.



A Diathermy Wand for just £185

Plus carriage and VAT



with a **FREE** pack of five tips worth **£50**

MIL 7501
Diathermy Wand 5mm Ø Re-usable

Diathermy Tips are available in packs of 5

MIL 7512
Diathermy Tips - single use, right angled tip - for MIL 7501



MIL 7513
Diathermy Tips - single use, 'J' tip - for MIL 7501



MIL 7514
Diathermy Tips - single use, spatula tip - for MIL 7501



MIL 7515
Diathermy Tips - single use, ball tip - for MIL 7501



MIL 7516
Diathermy Tips - single use, needle tip - for MIL 7501



SIGH

Surgical Instrument Group Holdings Limited 89a Gloucester Road Croydon CR0 2DN
 Telephone 020 8683 1103 Lo Call 0845 310 1103 Fax 020 8683 1105 Email info@sighltd.com



www.sighltd.com